

Vari-Strip™ & Vari-Plate™

PCR Tube Strips in Frames

Vari-Strips® and Vari-Plates™

PCR tube strips in frames - excellent flexibility


The 4titude® Vari-Strip™ and Vari-Plate™ system offers total flexibility in plate usage. It allows the user to insert or remove strips of 8 tubes from a 96 well plate frame.

- Vari-Strips™ can be used as a stand alone product or inserted in one of two Vari-Plate™ frames for ease of handling
- Vari-Plate™ frames are available for use on Roche LightCycler® 480 and for universal use
- Pre-loaded frames are available
- Sealable using cap strips, adhesive seals or heat seals



Figure 16: Vari-Strips™ can be used as a stand alone product or with a frame for ease of handling.

Vari-Strip™ 8 Well PCR Tube Strips, Low Profile



Features

- Thin-walled polypropylene tubes for optimum PCR and qPCR results
- Hole in end tab of white strips to visualise alphanumeric grid reference on the Vari-Plate™


Vari-Strip™ 8 Well PCR Tube Strips, Low Profile

Low profile, 0.1 ml PP wells, working volume: <100 µl, total well capacity: 200 µl

Code	Details	Strips/Case
4ti-0753	<input type="checkbox"/> W white wells	120
4ti-0753/C	<input type="checkbox"/> C clear wells	120

Code	Details	Tube and Cap Strips/Case
4ti-0754	<input type="checkbox"/> W white wells	plus strips of 8 flat optical caps 120
4ti-0754/C	<input type="checkbox"/> C clear wells	plus strips of 8 flat optical caps 120


Vari-Plate™ Frames & Vari-Plates™



Vari-Plate™ 96 Well Skirted Frame

PC, white, cut corner H1, universal use

Code	Details	Frames or Plates/Case
4ti-0757-F	Empty frame	10
4ti-0753/757	4ti-0757-F pre-loaded with Vari-Strips™, white	50
4ti-0753/C/757	4ti-0757-F pre-loaded with Vari-Strips™, clear	50



Vari-Plate™ 96 Well Semi-Skirted Frame, Roche Style

PC, white, cut corner H12, designed for use on Roche LightCycler® 480

Code	Details	Frames or Plates/Case
4ti-0950W-F	Empty frame	10
4ti-0753/950W	4ti-0950W-F pre-loaded with Vari-Strips™, white	50